



Serology Laboratory Bulletin January 19, 2011

Limited Quantity of Specimen from Potential Donor

When a specimen is submitted to NIT Serology Lab and is determined to be of such a volume that some tests may not be performed, the specimen is called "limited specimen". The "limited specimen" will be processed according to the following testing scheme unless the client directs otherwise. The testing will proceed in the sequence listed below.

First (1st) Priority Tests

1. Anti-HIV 1/2
2. HBsAg
3. Anti-HCV
4. Anti- HBC
5. Anti HTLV 1/II
6. Anti-HBs (if HBcore antibody is Positive)
7. ABO Rh

Second (2nd) Priority Tests

1. CMV IgG
2. CMV IgM
3. EBV (VCA) IgG
4. EBV(VCA) IgM
5. EBNA IgG
6. T cruzi Ab
7. Hgb A1C